



Inflammasome Signaling

Signal 1



Signal 2

Inflammasome activators are indicated in GREEN.




Key Products for Inflammasome Signaling

Visit www.cellsignal.com for full listing of antibodies, related kits, and reagents.


ANTIBODIES

	APPLICATIONS	REACTIVITY
#12948 AIM2 (D5X7K) Rabbit mAb	WB, IP	H
#67824 ASC (D2W8U) Rabbit mAb (Mouse Specific)	WB, IP, IF-IC	M
#3866 Caspase-1 (D7F10) Rabbit mAb	WB, IP	H
#4199 Cleaved Caspase-1 (Asp297) (D57A2) Rabbit mAb	WB, IP	H
#4450 Caspase-4 Antibody	WB	H
#46680 Caspase-5 (D3G4W) Rabbit mAb	WB, IP	H
#14340 Caspase-11 (17D9) Rat mAb	WB, IP	M
#14970 IF16 (D8B5T) Rabbit mAb	WB, IP, IHC-P	H
#12242 IL-1 β (3A6) Mouse mAb	WB, IHC-P	H, M
#12703 IL-1 β (D3U3E) Rabbit mAb	WB, IF-IC, F	H
#13627 IL-1 β (D3U3E) Rabbit mAb (PE Conjugate)	F	H
#31202 IL-1 β (D6D6T) Rabbit mAb (Mouse Specific)	WB, F	M
#83186 Cleaved-IL-1 β (Asp116) (D3A3Z) Rabbit mAb	WB, IP, IF-IC	H
#4990 NALP1 Antibody	WB	H, M, R, (Mk)
#15101 NLRP3 (D4D8T) Rabbit mAb	WB, IP	H, M, (Hm)
#12421 NLRC4 (D5Y8E) Rabbit mAb	WB, IP	H
#13833 TMS1 (E1E3I) Rabbit mAb	WB, IP	H

Untreated




LPS-treated




IL-1 β (D3U3E) Rabbit mAb #12703: Confocal IF analysis of THP-1 cells, untreated (left) or LPS-treated (500 ng/ml, 2 hr; right), using #12703 (green). Actin filaments were labeled with DyLight™ 554 Phalloidin #13054 (red).


Untreated




LPS-treated



Cleaved-IL-1 β (Asp116) (D3A3Z) Rabbit mAb #83186: Confocal IF analysis of THP-1 cells, differentiated with TPA #4174 (80 nM, 24 hr) and subsequently treated with (right) or without (left) Lipopolysaccharides (LPS) #14011 (1 μ g/ml, 6 hr), using Cleaved-IL-1 β (Asp116) (D3A3Z) Rabbit mAb (green). Actin filaments were labeled with DyLight™ 554 Phalloidin #13054 (red). Blue pseudocolor = DRAQ5® #4084 (fluorescent DNA dye).




Cleaved-IL-1 β (Asp116) (D3A3Z) Rabbit mAb #83186: WB analysis of extracts from cells or media collected from THP-1 cells, differentiated with TPA #4174 (80 nM, overnight) and subsequently treated with (+) or without (-) Lipopolysaccharides (LPS) #14011 (1 μ g/ml, 6 hr), using #83186.




NLRP3 (D4D8T) Rabbit mAb #15101: WB analysis of extracts from mouse bone marrow-derived dendritic cells (BMDC) and various cell lines using #15101 (upper) and β -Actin (D6A8) Rabbit mAb #8457 (lower).


BMDM Untreated




BMDM + LPS




BMDM + LPS/ATP



J774A.1



Raw 246.7



ASC (D2W8U) Rabbit mAb (Mouse Specific) #67824: Confocal IF analysis of mouse primary bone marrow-derived macrophages (BMDM), either untreated (A) or treated with LPS (50 ng/ml, 4 hr, B) or LPS followed by ATP (5 mM, 45 min, C), and J774A.1 (D) or Raw 246.7 (E) cells, using ASC (D2W8U) Rabbit mAb (Mouse Specific) (green). Blue pseudocolor = DRAQ5® #4084 (fluorescent DNA dye). Note the translocation of ASC to inflammasomes following stimulation with LPS and ATP (white arrows).

APPLICATIONS & REACTIVITY KEY: WB Western Blotting / IP Immunoprecipitation / IHC Immunohistochemistry / IF Immunofluorescence / F Flow Cytometry / ChIP Chromatin Immunoprecipitation / -IC Immunocytochemistry / -P Paraffin / -F Frozen / E-P Peptide ELISA / H human / M mouse / R rat / Hm hamster / Mk monkey / C chicken / Mi mink / Dm D. melanogaster / X Xenopus / Z zebra fish / B bovine / Dg dog / Pg pig / Sc S. cerevisiae / All all species expected / () 100% sequence homology

© 2016 Cell Signaling Technology, Inc. Cell Signaling Technology, and CST are trademarks of Cell Signaling Technology, Inc. U.S. Patent No. 5,675,063. All other trademarks are the property of their respective owners.